PL

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Warnings and safety information

- This product is intended solely for its intended use, i.e. to assist the user's mobility in accordance with the manufacturer's instructions. Do not use for medical, rehabilitation or therapeutic purposes unless the device has the appropriate certifications.
- Before first use, please read the user manual and all warnings provided by the manufacturer carefully.
- The exoskeleton may affect the user's balance and gait – take particular care when walking, climbing stairs and moving on uneven surfaces.
- Do not use the device if there is mechanical damage, wear and tear of structural components, or malfunction of the assistive systems.
- The product is not intended for use by children or persons with reduced ability to assess risk, unless they are under the constant supervision of a responsible person.
- Do not exceed the maximum permissible operating parameters specified by the manufacturer (e.g. maximum user weight, range of motion, operating time).
- When using the device, appropriate personal protective equipment must be worn, particularly in high-risk work environments.
- Do not make any unauthorised modifications to the device's structure or software – this may compromise safety and invalidate the manufacturer's liability.
- If you experience pain, discomfort, dizziness or other worrying symptoms, stop using the device immediately and consult a specialist.
- The device may contain electrical or electronic components – avoid contact with water unless otherwise specified by the manufacturer.

Information on correct use

- Before each use, check the technical condition of the device, including fastenings, straps and support systems.
- The exoskeleton should be properly fitted to the user's body in accordance with the manufacturer's instructions – incorrect fitting may lead to injury.
- The device should only be used under the conditions specified by the manufacturer (e.g. temperature, humidity, operating environment).
- The product must be regularly maintained in accordance with the manufacturer's recommendations, including cleaning and inspecting moving parts.
- For battery-powered models, use only original or manufacturer-recommended power sources.
- Do not leave the device switched on and unattended.
- Store in a dry place, secured against unauthorised access.

Additional precautions

- The product is not personal protective equipment unless it is explicitly marked as such and complies with the relevant standards.
- Do not use in potentially explosive atmospheres or in the presence of flammable substances unless the manufacturer permits such use.
- Transport the device in a manner that protects it from mechanical damage.
- All repairs must be carried out exclusively by an authorised service centre.
- Dispose of the end-of-life product in accordance with local regulations on waste electrical and electronic equipment (where applicable).
- The manufacturer and distributor shall not be liable for any damage resulting from use contrary to the instructions or the intended purpose of the product.

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